



Chilltherm[®] National

Fully Qualified +2°C to +8°C Insulated Shipping Solutions
for 24 to 36 hours

- +2°C to +8°C protection for up to 36 hours
- Qualification data available
- Tested to ambient temperatures of between -2°C to +25°C
- Common componentry across a number of sizes
- Selection of sizes available world wide
- Easy to dispose/recycle
- On-site training available


Chilltherm® National

- The Chilltherm® National qualified shippers form the short to moderate duration range of our National, Regional and Global temperature controlled packaging solutions.
- The range is manufactured at our advanced production centres worldwide and draws on our extensive experience in the field of cold chain shipping. Chilltherm® National has been designed to meet exacting specifications and has been fully qualified, with all components being manufactured to our strict ISO9001:2000 quality control systems.

- All of the shippers in the range feature a common bill of materials, with different configurations for hot and cold cycles, providing a single cost effective solution for both summer and winter shipments. Extensive supporting data can also be requested to demonstrate compliance with GDP.

For further information please contact us.

Chilltherm® temperature testing 24 hr road (National)



| | External dimensions (mm) | Volumetric weight (kg) | Available product dimensions (mm) | Available product volume (l) | Weight of system (kg) | Duration (hrs) |
|--------------------------|--------------------------|------------------------|-----------------------------------|------------------------------|-----------------------|----------------|
| Chilltherm® CTN2 | 351 x 266 x 241 | 3.8 | 235 x 64 x 132 | 2 | 1.9 | 24 > 36 |
| Chilltherm® CTN3 | 316 x 316 x 340 | 5.7 | 170 x 105 x 175 | 3 | 2.0 | 24 > 36 |
| Chilltherm® CTN11 | 470 x 470 x 410 | 15.1 | 280 x 205 x 200 | 11 | 6.1 | 24 > 36 |
| Chilltherm® CTN19 | 470 x 470 x 535 | 19.7 | 280 x 205 x 325 | 19 | 7.6 | 24 > 36 |
| Chilltherm® CTN23 | 571 x 515 x 441 | 21.6 | 365 x 250 x 250 | 23 | 9.5 | 24 > 36 |
| Chilltherm® CTN47 | 571 x 515 x 719 | 35.2 | 365 x 250 x 520 | 47 | 17.4 | 24 > 36 |

United Kingdom
 +44 (0) 1442 828664
enquiries@laminarmedica.co.uk
www.laminarmedica.com



Czech Republic
 +420 38 331 0000
enquiries@laminarmedica.cz
www.laminarmedica.cz



laminarmedica
 insulated shipping systems



Distributors

USA
 Saf-T-Pak Inc.
 +1 410-553-6066
sales@saftpak.com
www.saftpak.com



Canada
 Saf-T-Pak Inc.
 +1 780-486-0211
sales@saftpak.com
www.saftpak.com



Singapore
 Aeris Dynamics
 +65 6296 8368
sales@aeris.com.sg
www.aeris.com.sg

